Periodic Inspection protocol



Liko[™] overhead rail system

3EN111001 Rev. 13

System information	Customer reference	
Type of rail system:	Agreement No:	
S/N:	Name:	
Approved Max Load:	Address:	
Year of installation:	Post Code	

Liko[™] Overhead Rail System must be thoroughly inspected at least once per year. Inspection and service must be carried out by Hill-Rom authorized personnel. See "Installation condition checklist" on next page, to determine how to perform inspection points 2 and 5.

Periodic Inspections are available for: Liko™ Overhead Lifts (3EN191001), FreeSpan and FreeStand (3EN301001), Accessories (3EN601001) and Mobile Lifts (3EN371001) North America only: For non Liko™ Rail system, use "Non Liko® Rail Installation Periodic Inspection" (193725)

If printed: Make a colour print of this instruction.

INSTALLATION		Acti	Action required:			Not approved
1	Decal rail marking					
2	General inspection					
3*	Rails					
4	End stops					
5	Ceiling/Wall attachments					
ADDITIONAL COMPONENTS		N/A	S/N:	Action required:	Approved	Not approved
6*	Rail Switches					
7	Ultra Twist					
8	Transfer Motor Traverse					
9	Activity Rail					
10	Pool Side Lifter					
EN	VIRONMENTAL IMPACT					
11	Corrosive environments					
LOAD TESTING			Max Load applied:		Approved	Not approved
12 Maximum load rail system			Kg: Lbs:			
REC	REQUIRED MEASUREMENTS: N/A Dimension and unit		Approved	Not approved		
3* Joint Gap mm: Inch:		Inch:				
6B* Rail Switch play			mm: Inch:			
6C*	C* Turn Table clearance mm: Inch:					
Inspection sign off Complete inspection according to the instructions (page 2-11), fill in this page and sign below. Approval to use the overhead system Approved If the overhead system has one or more inspection points "Not approved", the system must not be used.						
Acti	Action required: Actions according to the inspection items "NOT APPROVED" should be performed immediately. After performed actions					

sign below. If anything is unclear or if you have questions, please contact Hill-Rom or your local Hill-Rom representative.

Contact information is to be found at www.hill-rom.com.

Inspection performed by:

Date:

Next inspection:

Inspection performed in accordance with ISO 10535:2006 Annex B- Periodic inspection

Enhancing outcomes for patients and their caregivers:



Installation condition checklist Liko™ overhead rail system


This checklist helps to determine ways to perform inspection points 2 and 5, depending on what environment the system is installed in.

All other inspection points has to be performed according to the instructions.

Environment:	Accessible fixing points: All attachments are visible	Partly accessible fixing points: <i>E.g. false ceiling with</i> <i>inspection hatches</i>	Non accessible fixing points: E.g. false ceiling
Non corrosive environment	 Inspect a minimum of 20% (at least two) of the fixing points. 	 Inspect a minimum of 20% (at least two) of the fixing points. If less than 20% can be inspected, extended maximum load test must be performed on non visible fixing points, according to section 12. 	 Extended maximum load test according to section 12.
Corrosive environment e.g. bathroom or other high humidity environment.	 Inspect all fixing points. Perform corrosion inspection according to section 11. 	 Inspect a minimum of 20% (at least two) of the fixing points. If less than 20% can be inspected, extended maximum load test must be performed on non visible fixing points, according to section 12. Perform corrosion inspection according to section 11. 	• Extended maximum load test according to section 12.
Chlorinated corrosive environment e.g. indoor pool	 Inspect all fixing points. Perform corrosion inspection according to section 11. 	 Inspect all fixing points. Perform corrosion inspection according to section 11. NOTE! Make inspection with inspection hatches or similar. 	 Inspect all fixing points. Perform corrosion inspection according to section 11. NOTE! Make inspection with inspection hatches or similar.


 Δ In this document, this warning symbol indicates that special care should be taken. If instructions are not followed there is a risk of serious injury.

Instructions for inspection points Liko[™] overhead rail system



INSTALLATION					
1 Decal rail marking					
 Inspect the S/N on the decal on the rail, if available Make sure the maximum load at the decal on the rail is equal to or higher than the maximum load for the lift unit installed in the system. Verify that there is a valid Installation Certificate for the overhead lift system. If no certificate is available a new needs to be created. 	In the class we way that the transformed of the class of				
2 General inspection					
 Inspect that the overhead lift system has no visible damage or deformations and that it is clean. Inspect that the distance between the installation points and the overhang match the instructions for the rail and the max load for installed lift unit. <i>Se Tables in Overhead System Installation Hand Book 3EN680001.</i> Use the appropriate inspection method for the bolts. (This is given by the supplier of the bolt.) 					
3 Rails					
 3A Primary rail (straight rail and traverse) Inspect that joint gap (A) is max 2 mm / 0,08 inch. Make sure that the rail joint is supported by a fixing or a proper joint section. 					
 3B Secondary Rail (traverse) Make sure the secondary rail (B) run freely, by placing a force (F) on one side. Listen for abnormal sounds from movable parts. Check traverse carrier bolts and nuts. 	(B) F				
4 End stops					
 Inspect that safety-through bolts with locking nuts are mounted. Inspect that the end stop is correctly mounted. 					


Liko[™] overhead rail system




6 Rail Switches

General Inspection Rail Switches

- Visually inspect to ensure that the switch are free from any sharp edges or deformities.
- Make sure bolts, nuts and joint screws are tightened.
- Make sure the switch is clean.
- Run the Side Rail Switch and Traverse Switch all the way in and out. Listen for peculiar noises and vibrations. Make sure the carriage runs easy when passing over the play between the rails.
- Check that the Turntable or Ultra Twist[™] system rotates freely on its bearings.
- Verify all functions on the handcontrol, or other options to operate the switch.
- Make sure cables, cords and mains power supply are intact without damages.
- Check connector pins on Multi station and connector rails on top of Multi-connector if applicable.





Instructions for inspection points Liko[™] overhead rail system

9 A	9 Activity Rail					
•	Verify presence of decal with model type and serial number.					
Ger •	neral Inspection Visually inspect to ensure that the Activity Rail is free from					
•	any sharp edges or deformities or damage. Move the Activity Rail all the way in the railsystem. Listen for peculiar noises and vibrations. Make sure the Activity Rail is clean.					
	·					
•	All bolts and safety clips on the Activity Rail is classified as safety critical and must be thoroughly checked for wear.					
10	Poolside Lifter					
For	Service and manintenance steel construction on Poolside Lifter-	contact Abus (www.abuscranes.com)				
Ger	neral Inspection					
•	Verify presence of decal with model type. Visually inspect to ensure that the Poolside Lifter is free from any sharp edges or deformities or damage. Move the Poolside Lifter all the way. Listen for peculiar noises and vibrations. Make sure the Poolside Lifter is clean.					
104	Poolside Lifter, for ground mounting					
•	Inspect the floor foundation Make sure all bots and nuts ar tightened. Inspect that the end stop is correctly mounted.					
	Maximum Load 200 kg / 440 lbs	al a state of the				
108	Poolside Lifter, for wall mounting	line_				
•	Inspect the wall foundation Make sure all bots and nuts ar tightened.					
•	Inspect that the end stop is correctly mounted.					
	Maximum Load 200 kg / 440 lbs					

Liko[™] overhead rail system

11 Environmental Impact – corrosive environments

Due to the environment an overhead system is installed in, components may be subject to corrosion. High temperature, high relative humidity, poor ventilation, presence of chlorine and different combinations of these factors, will affect the corrosion rate.

Depending on material type a corrosion attack can occur suddenly or in other cases form gradually. The corrosion rate and type of corrosion attack might be different in one area of the installation compared to another.

 Δ Fixing points classified as safety critical, installed in a corrosive environment such as indoor pool or bathroom, must be inspected. When a component has reached a certain stage of corrosion it might need to be replaced.

Note! Print out in color.

Check for visible severe corrosion and material loss and identify if components need to be replaced.

Galvanized steel

These pictures describes the evaluation method for all galvanized steel components.









1. New bolt

2. Acceptable


3. Replace




- 1. A galvanized steel component.
- 2. White rust on a component appears when the surface treatment corrodes.
- 3. Red rust appears when the actual steel has started to corrode. Corroding steel will result in material loss and should therefore be replaced.
- 4. A component covered in red rust is unfit for use.


Powder coated steel

These pictures describes the evaluation method for all powder coated steel components.



- 1. A powder coated steel component.
- 2. Local discoloration may occur in close proximity to corroding non-painted components. Stains on the painted surface is acceptable.
- 3. Cracks in the paint and red corrosion under the paint is a sign of corroding steel. Corroding steel will result in material loss and should therefore be replaced.
- 4. A component with peeling coating, bubbles in the paint and red corrosion under the paint is unfit to use.






Liko[™] overhead rail system

12 B Extended maximum load test

TRAVERSE SYSTEM

- Apply the maximum load for the overhead rail system.
- Place the carriage with load at the end stop of the secondary rail (1). Move the secondary rail, with a pause under each fixing point, from one end stop to the other end stop of the first primary rail.
- Continue by moving the load diagonally through the centre of the system over to the other side, as the dashed line shows.
- Now continue by moving the secondary rail with the applied load, from (2), with a pause under each fixing point, from one end stop to the other end stop of the second primary rail.
- Listen for unusual noises and vibrations. Check for any abnormalities such as deflection, movement or resistance.







Liko AB Nedre vägen 100 SE-975 92 Luleå, Sweden +46 (0)920 474700 www.hill-rom.com I www.liko.com



Enhancing outcomes for patients and their caregivers: