

# BAJ (MEDLINECARELINE®)



The battery pack BAJ has been specially developed for use with the JUMBO-system .

They are easy to exchange through an integrated snap system, and can easily be mounted on the mounting brackets

## Usage:

- Duty cycle: 10 % or 2 min . continuous use then 18 min . not in use
- Ambient temperatures: +5 °C to +40 °C
- Storage temperature: -15 °C to +40 °C
- BAJ is a part of the JUMBO system . It is compatible with CBJ1/CBJ2, CHJ2 and CBJ Care
- Flammability rating: UL94-V0
- Approvals: IEC60601-1, ANSI/AAMI ES60601, CAN/CSA-22 .2 No 60601-1

## Charging:

BAJ can be charged by

- Charger CHJ2
- Control box CBJ1/2, CBJ Care, COBO
- BAJ with integrated DC plug can also be charged by use of the external charger CH01

## Warning

Check at regular intervals that the ventilation hole is undamaged and intact . The construction of the ventilation stub permits battery gasses to get out, but it does not permit penetration of water



WE IMPROVE YOUR LIFE